

Select Agent Identification Chart


Select Agent	<i>Bacillus anthracis</i>	<i>Yersinia pestis</i>	<i>Francisella tularensis</i>	<i>Brucella</i> species	<i>Burkholderia mallei</i>	<i>Burkholderia pseudomallei</i>
Illness	Anthrax	Plague	Tularemia	Brucellosis, Malta fever	Glanders	Melioidosis
Signs and Symptoms	<u>Cutaneous</u> : eschar, edema <u>Gastrointestinal</u> : nausea, fever, abdominal pain, bloody diarrhea <u>Inhalational</u> : initial fever, malaise, some improvement, then abrupt respiratory distress	Fever, chills, headache, coughing up blood, buboes, sepsis	Ulceroglandular disease, sudden onset, chills, head cold, pharyngitis, cough, chest pain or tightness, possible skin lesions (cutaneous)	<u>Acute</u> : fever, sweats, headache, anorexia <u>Chronic</u> : mimic miliary TB	<u>Cutaneous</u> : nodules w/ localized lymphadenitis <u>Systemic</u> : broncho or lobar pneumonia, lesions in liver and spleen	<u>Acute</u> : pneumonia, high fever, dyspnea, chest pain <u>Subacute</u> : low grade fever, malaise, anorexia, weight loss <u>Chronic</u> : mimic miliary TB (Vietnamese time bomb)
Zoonotic	Yes	Yes	Yes	Yes	Yes	Yes
Transmission	Inhaling spores (e.g. from wool) or consuming undercooked meat from infected animals	Fleas, rodents, prairie dogs	Ticks, mosquitoes, fleas, contaminated hay, water, infected carcasses, infected animals, aerosolized particles, rabbits	Cattle (<i>B. abortus</i>), goats, sheep, camels (<i>B. melitensis</i>), pigs (<i>B. suis</i>), dogs (<i>B. canis</i>)	Equines (horses, mules, donkeys)	Soil, water (direct contact with), or aerosols from environmental sources
Growth Media	<small>Rec. Tests</small>	<small>Rec. Tests</small>	<small>Rec. Tests</small>	<small>Rec. Tests</small>	<small>Rec. Tests</small>	<small>Rec. Tests</small>
SBA	Growth 18 hr, ground glass	Growth 24-48 hr, gray-white to opaque	Scant to no growth 48-72 hr**	Growth 48-72 hr, white colonies	Growth 48 hr, gray, translucent colonies	Growth 24 hr, white colonies, wrinkled colonies > 72 hr
CHOC	Growth 18 hr	Growth 24-48 hr	Growth 48-72 hr, gray to grayish white colonies	Growth 48-72 hr	Growth 48 hr	Growth 24 hr
MAC	No growth	Growth 48 hr, clear, non-lactose fermenter	No growth	No growth	Poor to no growth 48 hr	Growth 24 hr
Gram stain, morphology, & size	* Positive Large rods, spores may be present 1-1.5 x 3-5 µm	* Negative Plump rods 1-2 x 0.5 µm	* Negative Coccobacillus, faintly staining 0.2-0.5 x 0.7-1.0 µm	* Negative, Coccobacillus, faintly staining 0.4 x 0.8 µm	* Negative Coccobacillus or small rod 1.5-3 x 0.5-1 µm	* Negative Straight or slightly curved rods 2-5 x 0.4-0.8 µm
Hemolysis	* Non-hemolytic	Non-hemolytic	Non-hemolytic	Non-hemolytic	* Non-hemolytic	* Non-hemolytic
Oxidase	+	-	-	+	Variable, most -	+
Catalase	+	+	+	+	+	+
Urease	-	-	-	+	variable	variable
Motility	-	-	-	-	-	+
Indole	-	-	-	-	-	-
Satellite	NR	NR	-	-	NR	NR
Growth @ 25° C	NR	+	NR	NR	NR	NR
Growth @ 42° C	NR	NR	NR	NR	-	+
Polymyxin B or Colistin	NR	NR	NR	NR	* resistant	* resistant
Amoxicillin clavulanate	NR	NR	NR	NR	* susceptible	* susceptible
Penicillin	NR	NR	NR	NR	* resistant	* resistant
β-lactamase	NR	-	+	-	NR	NR
ID panel misidentification		<i>Shigella</i> , <i>Y. pseudotuberculosis</i> , <i>Acinetobacter</i>	<i>Haemophilus influenzae</i> , <i>Aggregatibacter</i>	<i>Roseomonas gilardii</i> , <i>Kodcuria</i> spp.	<i>Burkholderia</i> spp., <i>Ralstonia pickettii</i>	<i>Chromobacterium violaceum</i> , <i>B. cepacia</i>
Occurrence	Endemic in the USA	Endemic in Idaho	Endemic in Idaho	Endemic in Idaho	Not endemic in USA	Not endemic in USA

NR: Not a recommended test for identification of specific select agent

**Initial culture may be positive on SBA; subsequent passages negative for growth on SBA

Culture Growth, Gram Stain, and Biochemical Tests

Bacillus anthracis



B. anthracis colonies at 24 h on SBA

Comma-shaped colonies





Tenacity on SBA



Gram stain

Yersinia pestis




Yersinia pestis growth on SBA at (A) 48 h, (B) 72 h, (C) 96 h, (D) 96 h "fried egg"


Giemsa stain (100X from blood culture): note *Y. pestis* bipolar appearance




Gram stain: note that bipolar staining may be poor




Francisella tularensis



Growth on SBA (A) 48 h, (B) 72 h



Growth on CHOC (A) 48 h, (B) 72 h



Gram stain

Brucella species


Growth on SBA at 48 h



Growth on CHOC at 48 h




Gram stain



Burkholderia mallei

Growth on SBA at 48 h



Gram stain



Motility




Non-Motile

Burkholderia pseudomallei


Growth on SBA at 24




Growth on SBA at 48 h



Gram stain



Oxidase



Positive


Motility





Motile

Bacillus anthracis Rule-Out Flowchart

Major Characteristics of *Bacillus anthracis*
Gram stain morphology: Large, Gram-positive rods
Spores may be found in cultures grown in 5% CO₂ or upon extended incubation but not usually in clinical specimens.
Colony morphology: Ground glass appearance, no hemolysis or pigment on SBA, no growth on MAC



Gram-positive bacillus



B. anthracis on SBA

Beta-hemolytic (i.e., exhibits complete clearing of blood cells on SBA)?

No

Catalase positive?

No


Bacillus anthracis is ruled out. Continue with routine identification.

Yes


Motile?

No

Bacillus anthracis not ruled out.
 Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.
Report: Possible *Bacillus anthracis* submitted to LRN Reference Laboratory.



B. anthracis is catalase positive.



B. anthracis is non-motile.

Yersinia pestis Rule-Out Flowchart



24 hours on SBA, 35°C

Major Characteristics of *Yersinia pestis*
Gram stain morphology: Gram-negative, plump rods, 0.5 x 1-2 µm
Colony morphology: Slow growing, pinpoint (1-2 mm), gray-white to opaque colonies on SBA after 24 hours; non-lactose fermenter on MAC/EMB; grows well at 25-28°C and 35-37°C.
Specimen is blood, sputum, or lymph node aspirate.



Gram-negative bacillus



48 hours on SBA, 35°C



72 hours on SBA, 35°C

Oxidase negative?
 Catalase positive?
 Indole negative?
 Urease negative?

No → *Yersinia pestis* is ruled out. Continue with routine identification.

Yes to all


Yersinia pestis not ruled out.
 Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.
Report: Possible *Yersinia pestis* submitted to LRN Reference Laboratory.



Y. pestis is oxidase negative.



Y. pestis is urease negative.



Y. pestis is catalase positive.

Warning: Some of the automated identification systems do not identify *Y. pestis* adequately. *Y. pestis* has been falsely identified as *Y. pseudotuberculosis*, *Shigella*, H₂S-negative *Salmonella*, *Acinetobacter*, and *Pseudomonas* spp.


Differentiation of other important *Yersinia* species

<i>Yersinia</i> species	Oxidase	Catalase	Urea	Indole
<i>Y. pseudotuberculosis</i>	Negative	Positive	Positive	Negative
<i>Y. enterocolitica</i>	Negative	Positive	Positive	Variable
<i>Y. frederiksenii</i>	Negative	Positive	Positive	Positive
<i>Y. kristensenii</i>	Negative	Positive	Positive	Variable
<i>Y. ruckeri</i>	Negative	Positive	Negative	Negative
<i>Y. pestis</i>	Negative	Positive	Negative	Negative




Francisella tularensis Rule-Out Flowchart

Major Characteristics of *Francisella tularensis*
Gram stain morphology: Aerobic, pleomorphic, minute (0.2 to 0.5 x 0.7 to 1.0 μm), faintly staining, Gram-negative coccobacillus
Colony morphology: No growth on MAC, scant to no growth on SBA after >48 hours
 Produces 1-2 mm gray to grayish-white colonies on CHOC after >48 hours.
Safety: Perform all work in a Class II biosafety cabinet using BSL-3 precautions.



Gram-negative coccobacillus




24 hours on SBA, 35°C




48 hours on SBA, 35°C



24 hours on CHOC, 35°C




48 hours on CHOC, 35°C



72 hours on CHOC, 35°C



F. tularensis is oxidase negative.



F. tularensis is weakly positive or negative on the catalase test.

Oxidase negative and catalase weakly positive or negative?

Yes

β-lactamase positive?
No growth on MAC?

Yes to both

Does the organism satellite around *S. aureus* on SBA?

No

No

Francisella tularensis is ruled out.
Continue with routine identification.

Yes

Francisella tularensis not ruled out.
Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.
Report: Possible *Francisella tularensis* submitted to LRN Reference Laboratory.

Warning: Automated identification systems may key out as non-*F. tularensis* (e.g., *Haemophilus influenzae* and *Aggregatibacter* spp.).

Brucella spp. Rule-Out Flowchart

Major Characteristics of *Brucella* species

Gram stain morphology: Small (0.4 x 0.8 µm), Gram-negative coccobacillus
THINK BRUCELLA

Colony morphology: Poor growth after 24 hours incubation on SBA and CHOC.
Incubate plates for at least 2 additional days if no growth in 24 hours.
Organism does **not** grow on MAC.


Safety: Perform ALL further work in a Class II biosafety cabinet using BSL-3 precautions when *Brucella* spp. is suspected.



Growth at 24 hours on SBA, 35°C



Growth at 48 hours on SBA, 35°C



Brucella spp. is oxidase positive.

Oxidase and catalase positive?

Yes to both

Urease positive?

Yes

Does the organism satellite around *S. aureus* on SBA?

No

Brucella spp. not ruled out.

Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.

Report: Possible *Brucella* spp. submitted to LRN Reference Laboratory.

Antimicrobial therapy: Rifampin or Streptomycin plus Doxycycline

No

Brucella spp. is ruled out. Continue with routine identification.

No

Yes



Gram-negative coccobacillus




Brucella spp. is urease positive.

Differentiation of *Brucella* from other urease positive, oxidase positive, Gram-negative coccobacilli

	<i>Brucella</i> spp.	<i>Psychrobacter immobilis</i>	<i>Paracoccus yeii</i>	<i>Psychrobacter phenylpyruvicus</i>	<i>Methylobacterium</i> spp.	<i>Oligella urolityca</i>	<i>Bordetella bronchiseptica</i> , <i>B. hinzii</i> , <i>Cupriavidus pauculus</i>	<i>Haemophilus</i> spp.
Gram stain morphology	Tiny coccobacillus, stains faintly	Coccobacillus/bacillus	Cocci in packets	Coccobacillus/bacillus, retains crystal violet	Vacuolated bacilli	Tiny coccobacillus	Bacillus	Tiny coccobacillus
Catalase	+	+	+	+	+	+	+	v
Oxidase	+	+	+	+	+	+	+	v
Urease	+	v	+	+	v	+	+	v
Motility	-	-	-	-	-	+(delayed)	+	-
SBA distinctions	-	Prefers 20°C, odor of roses	Mucoid	-	Pink, mucoid	-	-	No growth
MAC at 48 hr	-	-	-	-	-	-	+	-

Burkholderia mallei Rule-Out Flowchart




Growth at 24 hours on SBA, 35°C



Growth at 48 hours on SBA, 35°C




Growth at 72 hours on SBA, 35°C



B. mallei is catalase positive.

Major Characteristics of *Burkholderia mallei*
Gram stain morphology: Gram-negative coccobacilli or small rod
Colony morphology: Poor growth at 24 hours; better growth of gray, translucent colonies without pigment or hemolysis at 48 hours on SBA; poor or no growth on MAC at 48 hours; no distinctive odor



Gram-negative coccobacillus

Indole negative, catalase positive, non-hemolytic, poor growth or no growth on MAC?

Yes

Polymyxin B or colistin resistant (no zone)
 Amoxicillin-clavulanate susceptible
 Penicillin resistant

Yes

Scant or no growth at 42°C

Yes

Non-motile?

Yes

Burkholderia mallei is not ruled out.
 Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.
Report: Possible *Burkholderia mallei* submitted to LRN Reference Laboratory.
Additional screening test: *B. mallei* is arginine positive, unlike many other *Burkholderia* spp.


No

No

No

Burkholderia mallei is ruled out. Continue with routine identification.

No




B. mallei is non-motile.

Burkholderia pseudomallei Rule-Out Flowchart

Major Characteristics of *Burkholderia pseudomallei*

Gram stain morphology: Gram-negative rod, straight or slightly curved; may demonstrate bipolar morphology at 24 hours and peripheral staining, like endospores, when cultures are older

Colony morphology: Moderate growth at 24 hours, good growth of white colonies at 48 hours on SBA; may develop wrinkled colonies in time; non-pigmented; non-hemolytic; often demonstrates strong characteristic musty, earthy odor; growth on MAC at 48 hours.



Gram-negative bacillus




Growth at 24 hours on SBA, 35°C



Growth at 48 hours on SBA, 35°C



Growth at 72 hours on SBA, 35°C




B. pseudomallei is catalase positive.

Growth on MAC?

No → Rule out other agents, such as *Burkholderia mallei*, *Brucella* spp. and *Francisella tularensis*.

Yes
Oxidase positive, catalase positive, and indole negative?



B. pseudomallei is oxidase positive.

Yes to all
Polymyxin B or colistin resistant (no zone) or Growth on *B. cepacia* selective agars

No → *Burkholderia pseudomallei* is ruled out. Continue with routine identification.

Yes to one
No hemolysis on SBA
White to creamy yellow colonies on SBA

No → *Burkholderia pseudomallei* is ruled out. Consider *Chromobacterium violaceum* or indole-negative *Vibrio* spp.

Yes

Burkholderia pseudomallei is not ruled out.
Call Idaho Bureau of Laboratories (208-334-2235) to notify LRN-B personnel and send suspected agent according to applicable shipping regulations.
Report: Possible *Burkholderia pseudomallei* submitted to LRN Reference Laboratory.
Additional screening test: *B. pseudomallei* is arginine positive, unlike other *Burkholderia* spp. Unlike *B. mallei*, *B. pseudomallei* grows at 42°C at 48 hours and is motile.