## TRIPLE PACKAGING FOR SPECIMEN SHIPMENT

Category A and Category B Infectious Shipments











## PRIMARY LEAKPROOF SPECIMEN CONTAINERS

Labeled with appropriate patient information.

### **BUBBLE WRAP OR OTHER CUSHIONING MATERIAL**

Wrap individual primary specimen containers to prevent contact between specimen containers.

# ABSORBENT SHEET BETWEEN PRIMARY SPECIMEN CONTAINER AND BIOHAZARD BAG

Amount of absorbent must be enough to absorb contents.

SECONDARY PACKAGING
LEAKPROOF BIOHAZARD BAG

TYVEK ENVELOPE, PRESSURE TESTED TO 95 kPa

#### ALL OF ABOVE GO INTO OUTER BOX FOR SHIPPING

Use cushioning material to reduce specimen movement within the outer container!