

MERCURY, LEAD AND CADMIUM IN WHOLE BLOOD ANALYSIS REQUEST FORM

COMPLETE ALL INFORMATION

Submitter Name:		Submitter Organization:	
Mailing Address:			
City:		State:	Zip:
Phone:		E-mail Address:	
Secure Fax Number:			

State of Idaho
Bureau of Laboratories
2220 Old Penitentiary Rd.
Boise, ID 83712
208-334-2235



Customer Specimen ID	Collected by	Collection Date	Collection Time	Gender	Date of Birth	Tobacco Use Frequency	Known Chem Exposure?	Symptomatic Chem Exposure?	Laboratory Use Only		
									Lab Sample #	Bottle Type	Preservative

Special Instructions: _____

General Information

Optimal amount of specimen is 1-2 mL, minimum is 0.25 mL. Specimens must be collected in evacuated blood tubes with EDTA anticoagulant (tubes with royal blue closures, certified for heavy metal/trace element analysis, are highly recommended). Draw the blood through a stainless steel needle. If more than one evacuated tube of blood is to be drawn from an individual, the trace metals tube should be drawn second or later. Submit two blank (empty) containers per lot of specimen containers. Handle and store blank containers identical to specimens for testing. Store and transport specimens at 5 ± 3 °C. If shipping specimen to IBL, package for transport according to IATA regulations.

LABORATORY USE ONLY

Temperature: _____ Shipper: _____ Number of Bottles / Sample: _____ Number of Samples / Order _____ Lab Order ID: _____

Date Received: _____ Time Received: _____ Received By _____ Revised 11/28/07