



IDAHO DEPARTMENT OF HEALTH & WELFARE
DIVISION OF PUBLIC HEALTH

Test Title	Lyme Disease Antibody Methodology: EIA with Reflex to Western Blot Confirmatory Testing
Sample Requirements	<ol style="list-style-type: none">1. Sample type(s): 2 ml serum and a completed Lyme disease report form (call virology at 334-2235 for a faxed copy)2. Rejection criteria: N/A
Sampling Materials	<ol style="list-style-type: none">1. Sample container: blood tube or serum transfer screw cap tube
Procedural Notes	<ol style="list-style-type: none">1. Clinical Test Request Form2. Referred to CDC, Fort Collins, Colorado
Shipping Instructions	<ol style="list-style-type: none">1. Temperature/preservative instructions: room temperature2. Package according to Biological Substance, Category B, shipping guidelines.3. Ship to: Idaho Bureau of Laboratories ATTENTION: Virology Laboratory 2220 Old Penitentiary Rd Boise, ID 83712
Reporting and Turnaround Time (TAT)	<ol style="list-style-type: none">1. TAT: 2-6 weeks2. This disease must be reported to your local health district or to the state Bureau of Communicable Disease Prevention according to the rules and regulations governing Idaho reportable diseases (IDAPA 16.02.10).3. Reference range: N/A