



IDAHO DEPARTMENT OF HEALTH & WELFARE
DIVISION OF PUBLIC HEALTH

Test Title	<i>Aeromonas</i> spp. Detection In Stool Methodology: Culture
Sample Requirements	<ol style="list-style-type: none">1. Sample type(s): stool in parapak C&S, Modified Carey-Blair or equivalent - fill to line (approximately 5 ml)2. Unpreserved stool samples must be received within 2 hours of collection
Sampling Materials	<ol style="list-style-type: none">1. Sample container: Stool parapak C&S container or screw cap specimen cup
Procedural Notes	<ol style="list-style-type: none">1. Clinical Test Request Form2. Restrictions: done for investigation of foodborne illness only. Contact your local health district before requesting this test3. CPT code: 87046
Shipping Instructions	<ol style="list-style-type: none">1. Temperature/preservative instructions: hold and ship preserved samples at room temperature. Transport fresh samples at refrigeration temperature.2. Package according to Biological Substance, Category B, shipping guidelines.3. Ship to: Idaho Bureau of Laboratories ATTENTION: Bacteriology Laboratory 2220 Old Penitentiary Rd Boise, ID 83712
Reporting and Turnaround Time (TAT)	<ol style="list-style-type: none">1. TAT: 6 business days2. Reference range: N/A