



Student Name _____
 Date _____
 Evaluator _____

Skill Sheet

Nerve Agent Antidote – EMR, EMT, AEMT

Skill		Pass	Needs Improvement (Please comment)
1. Takes BSI precautions, may include specialized PPE			
2. Gathers necessary equipment			
3. Assures medication is indicated for the patient			
4. Determines if patient has taken any doses			
5. Checks medication for expiration date			
6. Contacts medical direction for authorization, if administration is not covered in standing orders			
SPECIFIC MEDICATIONS			
Mark-1™	DuoDote™		
7. Explains procedure to patient	7. Explains procedure to patient		
8. Tears open plastic pouch at any of the notches and removes the Mark-1 Auto-Injector	8. Tears open plastic pouch at any of the notches and removes the DuoDote Auto-Injectors		
9. Places Mark-1 auto-injector 1 (atropine-smaller one) in dominant hand and grasps firmly with Black tip pointing downward	9. Places DuoDote in dominant hand a grasps firmly with Green tip pointing downward		
10. Removes gray safety cap from auto-injector 1	10. Pulls off gray safety release, taking care to not touch the green tip		
11. Keeps fingers clear of both ends of the auto-injector	11. Keeps fingers clear of both ends of the auto-injector		
12. Firmly pushes the black end of the auto-injector against the lateral side of the patient's thigh, midway between waist and knee	12. Finds injection site on mid-outer thigh area		
13. Continues to push firmly until they feel the auto-injector trigger	13. Firmly pushes Green Tip strait down (at 90° angle), continuing until they feel auto-injector trigger		
14. Holds the auto-injector in place until the medication is injected (abt 10 sec) Massages the site for several seconds	14. After the DuoDote Auto-injector triggers, firmly holds in place for 10 seconds. Massages the site for several seconds		
15. Removes the Mark-1 Auto injector from thigh and inspect tip for needle	15. Removes the DuoDote auto-injector from the thigh and inspects tip for needle		



Student Name _____
 Date _____
 Evaluator _____

Skill Sheet

Nerve Agent Antidote – EMR, EMT, AEMT

16. If needle is NOT visible, makes sure Gray Safety Release is removed and repeats injection procedure	16. If needle is NOT visible, makes sure the Gray Safety Release is removed and repeats injection process		
17. Pushes needle against a hard surface to bend needle back against auto-injector, then puts back in plastic pouch	17. Pushes exposed needle against hard surface until it bends back, then puts back in plastic pouch		
18. Candidate repeats steps 9-17 for auto-injector 2 (pralidoxime- large one)		/9	
19. Keeps used auto-injectors with patient	19. Keeps used auto-injector(s) with patient		
Total Steps Passed			
Minimum steps needed to pass (assuming no critical fails below)			

Critical Criteria (Automatic Fail)

- _____ Did not take BSI precautions
- _____ Did not determine medication was indicated for patient
- _____ Administered antidote in any way that is harmful to patient or provider
- _____ Was not able to administer the antidote
- _____ Did not take care of needle properly

Final Score (circle one)	Pass	Fail
---------------------------------	-------------	-------------